Page 1 of 1

CA-8-65-28861-00166

Certificate of Calibration

Detail: BLOOD PRESSURE ID No.: BPA-NS-02

Manufacturer: ALRK2 Model: STANDBY Serial No.: -

Department: โรงพยาบาลหนองหิน Section: รพ.สต.น้อยสามัคคื

หนองหิน จังหวัดเลย

Room Ambient: Temperature: 23 °C (15-25 °C) Humidity: 50 %RH (20-85 %RH)

Measurement Method: Follow the blood pressure monitor oscillometric calibration method of Wo-CAL-BP-co. By using

measuring instruments that have been calibrated by a standardized agency and calibrated

Cert. No.:

under specified environmental conditions.

Standard Used:

Equipment	Manufacture	Model	Serial No.	Certification No.	Cal Date
Pressure Calibrator	Fluke	717 100G	3366021	21P2453	7/15/2021

Measurement Results:

UUC Setting (mmHg)	STD Reading (mmHg)	Error (mmHg)	Uncertainty (+/-mmHg)	Accept value (+/-mmHg)
50.00	50.00	0.00	2.33	4.00
100.00	100.00	0.00	2.33	4.00
120.00	120.00	0.00	2.33	4.00
150.00	150.00	0.00	2.33	4.00
180.00	180.00	0.00	2.33	4.00
200.00	200.00	0.00	2.33	4.00

The reported uncertainty is based on a standard uncertainty multiplied by a coverage factor

k=2, Providing a level of confidence of approximately 95%

Cal by : Mr.Kitipong Saengkaew Approved by :

Plut. T.

(Mr.Phatchara Temrat)

Biomedical Engineer

Comment: This report certifies the calibrate results as of the date, the location and the conditions of calibrate only.